ACCU-CHEK, ACCU-CHEK TENDERLINK and AKKY-CHEK are trademarks of Roche.

© 2021 Roche Diabetes Care

Made in Mexico

Unomedical a/s  
Aaholmvej 1-3, Osted  
4320 Lejre, Denmark  
[www.infusion-set.com](http://www.infusion-set.com)

Last update: 2020-06


**Distributed by:**  
Roche Diabetes Care GmbH  
Sandhofer Strasse 116  
68305 Mannheim, Germany

03692051001/R

**ACCU-CHEK® TenderLink**

**Infusion set**  
Instructions for Use

**ACCU-CHEK®**

Instructions for use	4	Návod na použitie	230
Gebrauchsanweisung	20	הוראות שימוש	261
Instructions d'utilisation	36	Használati útmutató	262
Istruzioni per l'uso	52	Navodila za uporabo	278
Instrucciones de uso	68	Instrucțiuni de utilizare	294
Instruções de utilização	84	Инструкции по использованию	310
Gebruiksaanwijzing	100	Інструкція із застосування	326
Brugsanvisning	116	Kasutusjuhend	342
Bruksanvisning	132	Lietošanas pamācība	358
Bruksanvisning	148	Naudojimo taisyklės	374
Käyttöohje	164	إرشادات الاستخدام	405
Οδηγίες χρήσης	180		
Instrukcja obsługi	196		
Návod k použití	212		

<b>Instructions for use</b>	<b>en</b>
<b>Gebrauchsanweisung</b>	<b>de</b>
<b>Instructions d'utilisation</b>	<b>fr</b>
<b>Istruzioni per l'uso</b>	<b>it</b>
<b>Instrucciones de uso</b>	<b>es</b>
<b>Instruções de utilização</b>	<b>pt</b>
<b>Gebruiksaanwijzing</b>	<b>nl</b>
<b>Brugsanvisning</b>	<b>da</b>
<b>Bruksanvisning</b>	<b>no</b>
<b>Bruksanvisning</b>	<b>sv</b>
<b>Käyttöohje</b>	<b>fi</b>
<b>Οδηγίες χρήσης</b>	<b>el</b>
<b>Instrukcja obsługi</b>	<b>pl</b>
<b>Návod k použití</b>	<b>cs</b>

## ***ACCU-CHEK® TenderLink***

### Infusion set

#### **Indications for use**

The Accu-Chek TenderLink infusion set is indicated for subcutaneous infusion of insulin administered by an external pump. The infusion set is indicated for single use.

#### **Intended use**

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes. The infusion set has been tested to be compatible with insulin approved for subcutaneous infusion.

#### **Product description**

The Accu-Chek TenderLink infusion set consists of a cannula and tubing, i.e. cannula (head set, picture 1a) and tubing (transfer set, picture 1b).

## Contraindications

Accu-Chek TenderLink infusion set is neither intended nor indicated for use with blood or blood products (I.V. infusion).

### WARNINGS

- Carefully examine the instructions for use before inserting Accu-Chek TenderLink infusion set. Failure to follow instructions may result in pain or injury.
- The Accu-Chek TenderLink infusion set is only sterile and non-pyrogenic if the packaging is undamaged and unopened. Do not use it if the packaging is already open or is damaged.
- When using Accu-Chek TenderLink infusion set for the first time, do so in the presence of a healthcare professional.
- Inaccurate medication delivery, site irritation or infection may result from improper insertion and maintenance of infusion site.
- Do not re-insert the introducer needle into the infusion set. This could cause a tearing or damage of the soft cannula and may lead to unpredictable medication flow.

- Reuse of the infusion set may cause site irritation or infection and may lead to inaccurate medication delivery.
- Be sure that the needle cover is removed before insertion is made.

### **PRECAUTIONS**

- Never try to fill or free a clogged tubing while the infusion set is inserted. This may result in unpredictable medication flow.
- Do not leave air in the infusion set. Make sure to fill the tubing completely. Consult pump manufacturer's instructions for use.
- Accu-Chek TenderLink infusion set is a single use item, which must be disposed of after usage. Do not clean or re-sterilise.
- Check the infusion set frequently to ensure that the soft cannula remains firmly in place. Replace with a new infusion set at a new site, if the infusion set is not in place. Since the cannula is soft, it will not cause any pain if it slips out, and this may happen without notice. The soft cannula must always be completely inserted to receive the full amount of medication.
- Do not in any way put disinfectants, perfumes, deodorants or other products containing disinfectants in contact with the luer lock connector or the tubing. This may impact the integrity of the infusion set.

- Replace the infusion set if the adhesive pad becomes loose.
- Replace the infusion set every two to three days, or per your healthcare professional's instructions.
- If the infusion site becomes irritated or infected, replace the infusion set and use a new site until the first site has healed.
- Use aseptic techniques when temporarily disconnecting Accu-Chek TenderLink infusion set. Consult with healthcare professional on how to compensate for missed medication when disconnected.
- People with small or large amounts of subcutaneous fat should choose the insertion angle carefully, as the cannula could be placed in underlying muscle or dermal layer which may lead to unpredictable medication delivery. Normal insertion angle is between 20–45°. Consult your healthcare professional concerning this matter.

### **Recommendations**

- Check your blood glucose level one to three hours after insertion of the Accu-Chek TenderLink infusion set. Measure your blood glucose on a regular basis. Discuss this with your healthcare professional.
- Do not change infusion set just prior to bedtime, unless blood glucose can be checked one to three hours after insertion.



- If your blood glucose level becomes unexplainably high, or occlusion alarm occurs, check for clogs and leaks. If in doubt, change infusion set since soft cannula could be dislodged. Discuss plan to bring down your blood glucose level with your healthcare professional should any of these problems arise. Test blood glucose level to make sure problem is corrected.
- Carefully monitor blood glucose levels when disconnected and after reconnecting.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

### **Storage & disposal**

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Dispose of the introducer needle in a proper sharp container according to local laws, to avoid risk of needle stick.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

en


9


## Preparation and insertion

### Contents

- A Needle cover
- B Soft cannula
- C Adhesive pad
- D Introducer needle
- E Cannula casing
- F Tubing
- G Connector with connector needle
- H Circular protective cap
- I Disconnect cover


Wash your hands.


Suggested insertion sites. Choose an insertion site as recommended by your healthcare professional.

en



Clean the area of insertion with a disinfectant as directed by your healthcare professional. Leave skin to air dry for at least 30 seconds.


Remove the front paper from the adhesive pad.

6 Gently twist and pull to remove the needle cover. Make sure the soft cannula does not exceed the introducer needle.


7 Pinch the skin and insert the infusion set at a 20–45-degree angle.


8 Place a finger in front of the cannula casing and remove the introducer needle by pressing on the grips on the side of the introducer needle and pulling straight out.

en


9  
Keep your finger in front of the cannula casing and remove the back paper from the adhesive pad. Gently massage the adhesive pad onto the skin.


10  
Remove paper from the tubing and circular protective cap from the tubing connector needle. Fill the tubing until medication dispenses from the tip.

11

Consult pump manufacturer's instructions for use when connecting infusion set to the insulin pump.

12  
Place a finger in front of the cannula casing and push the connector needle straight in.

13  
Click.


14

Fill the empty space in the cannula with the following fill volume:  
U100 insulin:  
0.7 U (0.007 mL)


en

15


## Disconnect

Place a finger in front of the cannula casing. Disconnect the connector needle by pressing the grips on the side and pulling straight out.


Connect the circular protective cap to the connector needle.


Connect the disconnect cover to the cannula casing.


## Reconnect

Remove circular protective cap from the connector needle and if necessary, fill the tubing.

Remove disconnect cover from the cannula casing. Place a finger in front of the cannula casing and push the connector needle straight in.

Click.

en

Distributed in the United Kingdom by:  
Roche Diabetes Care Limited  
Charles Avenue, Burgess Hill  
West Sussex, RH15 9RY, **United Kingdom**

Accu-Chek Pump Careline <sup>1)</sup>:

UK Freephone number: 0800 731 22 91

ROI Freephone number: 1 800 88 23 51

<sup>1)</sup> calls may be recorded for training purposes

Some mobile operators may charge for calls to these numbers.

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